


Power Settings Alerion Main Computer


The computer that is holding your Alerion data needs to have its power options altered so the computer does not go to sleep or turn off its network card.

Sleep Mode

- Click Start
- Click Control Panel
- Click Power Options
- Click Change Plan Settings next to the plan that has the radio button lit.




- Click the drop down box next to 'Put the computer to sleep' and change it to Never.




- Click Save changes

Network Card

- Click Start
- Click Control Panel
- Click System
- Click Device Manager
- Click Network Adapters
- Right- Click on your Network card then choose properties on the little window that appears.



- Look for a Power Management tab on this windows. If you don't have one you are done.



- Click the Power Management tab and turn off 'Allow the computer to turn off this device to save power'.
- Click Ok and exit.